

# Emergency Disconnect Procedure

Model AK 96  
For Patients  
Graft or Fistula

1



**PRESS BLOOD PUMP BUTTON.**


Blood pump will stop and button will start flashing.

2




**USING THE PRODUCTS FOUND IN YOUR EMERGENCY KIT...**

3



**CLOSE TWO CLAMPS ON NEEDLES AND BLOODLINES. ADD SAFETY CLAMPS BELOW AND ABOVE THE CONNECTION.**

4



**DISCONNECT LINES AT THE CONNECTION SITE BY TWISTING.**

5



**GO TO A SAFE PLACE WITH YOUR EMERGENCY KIT.**

Produced by

**BCRenal**

Provincial Health Services Authority

