


# Emergency Disconnect Procedure

Model AK 98  
For Patients  
Catheter

**1**




**PRESS BLOOD PUMP BUTTON.**  
Blood pump will stop and button will start flashing.

**2**




**USING THE PRODUCTS FOUND IN YOUR EMERGENCY KIT...**

**3**



**CLOSE CLAMPS ON THE CATHETER AND BLOODLINES. ADD SAFETY CLAMPS BELOW AND ABOVE THE CONNECTION.**

**4**



**DISCONNECT LINES AT THE CONNECTION SITE BY TWISTING.**

**5**



**GO TO A SAFE PLACE WITH YOUR EMERGENCY KIT.**

